

## Model 3276HFM Low Cost Hand & Foot Monitor

### Features

- Measure Alpha, Beta, and Gamma Contamination
- Simple Setup and Use
- User-Adjustable Alarms
- User-Programmable Units of Measure, Either U.S. or SI Units
- Power: Standard Alkaline (4) "AA" Batteries or 9 Vdc Wall Mount Transformer

### Introduction

The Model 3276 is a versatile instrument that can be used for multiple radiation detection purposes. In this configuration, the Model 3276 is connected to both a Model 44-25 hand monitor detector and a Model 44-26 foot monitor detector to monitor personnel for alpha, beta, and gamma contamination. An optional Model 44-9 detector can also be included for frisking.

### Configurations

<b>Product</b>	 <b>Model 3276HFM</b>	 <b>Model 3276HFM with Frisker</b>	 <b>Model 3276HFM with Frisker</b>
<b>Part Number</b>	48-4160-2	48-4160-6	48-4160-1
<b>3276 Meter</b>	✓	✓	✓
<b>44-25 Hand Detector</b>	✓	✓	✓
<b>44-26 Foot Detector</b>	✓	✓	✓
<b>Mount</b>	Stand	Stand	Wall
<b>44-9 Frisker</b>		✓	✓
<b>Detector Switch</b>	View hand counts, foot counts, or hand + foot counts	View hand counts, foot counts, or frisker counts	View hand counts, foot counts, or frisker counts
<b>Size (L x W x H)</b>	32 x 30.5 x 138 cm (12.6 x 12 x 54.3 in.)	32 x 30.5 x 138 cm (12.6 x 12 x 54.3 in.)	See Component Specifications

Ludlum Measurements, Inc. P.O. Box 810, Sweetwater, Texas 79556

Web: <https://www.ludlums.com> Tel: 800-622-0828 / 325-235-5494 / Fax: 325-235-4672 / Email: [sales@ludlums.com](mailto:sales@ludlums.com)

Note: specifications subject to change without notification. We are not responsible for errors or omissions.

Mar 2024

## Model 3276 Multipurpose Digital Meter

**Part Number:** 48-4160

Features include user-adjustable alarms and user-programmable units of measurement in either U.S. or SI units.

**DISPLAY:** 3-digit LCD with 20 mm (0.8 in.) digits

### DISPLAY RANGE:

- 0.00 cps to 999 kcps
- 0.00 cpm to 999 kcpm
- 0.00 Bq to 999 kBq
- 0.00 dpm to 999 kdpm
- 0.00  $\mu$ R/h to 999 R/h
- 0.00  $\mu$ Sv/h to 999 Sv/h

NOTE: Display Ranges should be set up by qualified personnel according to detector specifications

### USER CONTROLS:

- ON/OFF – Extra-long press to turn ON; Short press to acknowledge alarms; Extra-long press plus 3 seconds to turn OFF
- UNITS – Short press changes the units from count rate, to dose/exposure, or disintegration, extra-long press to change count time in COUNT Mode
- AUDIO – Short press to adjust audio level
- MODE – Short press alternates between modes: RATE (displays count rate), MAX (captures peak rate), and COUNT (captures rate or dose readings for a user-selectable preset count time from 0 to 10 minutes)

**RESPONSE TIME:** User-selectable from 1 to 60 seconds, or Auto-Response Rate FAST or SLOW

**ALARMS:** Two adjustable count rate, exposure/dose, and scaler alarm setpoints adjustable over the display range

**LOSS OF COUNT PROTECTION:** After a user-settable time interval (default 60 seconds) of no pulses from detector, the unit will flash a zero reading and the alarm audio will be triggered

**AUDIO:** Approximately 4.5 kHz, click audio greater than 65 dB at 0.6 m (2 ft), alarm audio greater than 72 dB

**POWER:** 4 alkaline "AA" batteries with approximately 50 hours expected life, or 9 Vdc input via wall mount transformer (input 100-250 Vac, 50-60 Hz)

**BATTERY LIFE:** Approximately 50 hours with the backlight ON, 16-hour low-battery warning

**ENVIRONMENTAL RATING:** NEMA rating of 2, IP rating of 51

**SIZE (H x W x L):** 5.8 x 13.5 x 19.1 cm (2.3 x 5.3 x 7.5 in.)

**WEIGHT:** 0.9 kg (2 lb) (excluding detectors and clips)

## Model 44-25 Hand Detector

**Part Number:** 47-1508

**APPLICATION:** Alpha, beta, gamma survey of hands

**DETECTOR:** 2 pancake-type, halogen-quenched GM tubes

**WINDOW (ACTIVE):** 31 cm<sup>2</sup> (4.8 in<sup>2</sup>)

**SENSITIVITY:** typically 6600 cpm per mR/hr

**SIZE (H x W x L):** 24.1 x 15.2 x 19.8 cm (9.5 x 6 x 7.8 in.)

## Model 44-26 Foot Detector

**Part Number:** 47-1509

**APPLICATION:** Alpha, beta, gamma survey of foot/shoe

**DETECTOR:** 3 pancake-type, halogen-quenched GM tubes

**WINDOW (ACTIVE):** 46 cm<sup>2</sup> (7.1 in<sup>2</sup>)

**SENSITIVITY:** typically 9900 cpm per mR/hr

**SIZE (H x W x L):** 27.4 x 17.8 x 22.9 cm (10.8 x 7 x 9 in.)

## Model 44-9 Frisker (Optional)

**Part Number:** 47-1539

**APPLICATION:** Alpha, beta, gamma survey

**DETECTOR:** Pancake-type, halogen-quenched GM tube

**WINDOW (ACTIVE):** 15.5 cm<sup>2</sup> (2.4 in<sup>2</sup>)


**SENSITIVITY:** typically 3300 cpm per mR/hr

**SIZE (H x W x L):** 4.6 x 6.9 x 27.2 cm (1.8 x 2.7 x 10.7 in.)

## Options

**Stand Wheel Option** (PN 4224-025)

For Model 3276HFM with stand mount, adds wheels to the base to make it easier to move the monitor between different locations.



See our website for additional options for the Model 3276.