

Model L-989-415

PET L-Block Shield



Ludlum Measurements, Inc.



Model L-989-415 PET L-Block Shield

Part Number: 99-9305

The 3.8 cm (1.5 in.) thick lead L-Block Shield provides exceptional protection to the clinician loading iodine and iridium.

The face shield is made of optically clear lead glass, cantilevered for unimpaired viewing of the work area. The leaded glass has a density of 5.2 g/cm³ to give maximum eye and face protection.

LEAD GLASS: 2 cm lead equivalent (¹³⁷Cs energy)

LEAD GLASS DIMENSIONS: 20.3 x 30.5 x 5.1 cm 8 x 12 x 2 in (H x W x D)

LEAD DENSITY: 11.35 g/cm²

FRONT SHIELD DIMENSIONS: 28.3 x 31.5 cm (11.1 x 12.4 in.) (H x W)

FINISH: Durable tan textured

OVERALL DIMENSIONS: 42.9 x 31.5 x 34.3 cm (16.9 x 12.4 x 13.5 in.) (H x W x D)

WEIGHT: 95.3 kg (210 lb)