

Model L-929-10 and L-929-91

Syringe Holder/Pig and Syringe Carrier



Ludlum Measurements, Inc.



Part Number: 99-9301



Part Number 99-9300

Model L-929-10 Syringe Holder

The Model L-929-10 Syringe Holder was designed to accommodate unshielded syringes that contain isotopes ranging from 1 cc to 20 cc. Lead, 1.27 cm (0.5 in.) thick, surrounds the well and provides protection for the user. The extra-wide diameter base provides stability to prevent accidental tipping or damage from jostling.

OUTSIDE DIMENSIONS:

16.5 x 9.8 cm
(6.5 x 3.9 in.)
(H x W)

INSIDE DIMENSIONS:

14.3 X 2 cm)
(5.6 X 0.8 in.)
(H x W)

WEIGHT: 2.7 kg (6 lb)

Model L-929-91 Syringe Carrier

Syringe Carrier, Model L-929-91, is constructed of stainless steel with 3.2 mm (0.13 in.) lead completely enclosed between the inner and outer walls. It is designed to safely contain syringes, vials, and ampules up to 22.9 cm (9 in.) length. The device may be used for storage and transport. The lid opens 180° for accessibility, overlaps for security, and has a latch to prevent streaming.

OUTSIDE DIMENSIONS:

6.4 x 6.9 x 24.4 cm
(2.5 x 2.7 x 9.6 in.)
(H x W x L)

INSIDE DIMENSIONS:

5.1 X 5.1 X 22.9 cm
(2 x 2 x 9 in.)
(H x W x L)

WEIGHT: 3.6 kg (8 lb)