

# Model L-929-47 and Model L-912-00

## Lead Lined Storage/Waste Container and Lead Bricks



Ludlum Measurements, Inc.



### Model L-929-47 Lead Lined Storage/ Waste Container

Part Number: 99-9302

The Lead Lined Storage Container is ideal for storing syringes or other items that have been contaminated with low-level energy gamma residue. The 0.32 cm (0.13 in.) lead shielding allows contaminated items to be stored until they are properly decayed.

#### OVERALL DIMENSIONS:

16.5 x 12.7 cm  
(6.5 x 5 in.)  
(H x Dia)

#### SHIELDING:

0.3 cm  
(0.125 in.)

**WEIGHT:** 3.2 kg (7 lb)

### Model L-912-00 Lead Bricks

Part Number 99-9580

Ludlum Lead Bricks are ideal for building temporary or permanent shielding for hot labs or isotope storage/waste facilities. The brick ends are saw cut so that they will fit together tightly, eliminating open spaces between stacked bricks.

Lead Bricks are available in two sizes:

#### Standard:

5.1 x 10.2 x 20.3 cm  
(2 x 4 x 8 in.)  
(D x W x L)  
Part Number 99-9580

#### Optional:

5.1 x 10.2 x 15.2 cm  
(2 x 4 x 6 in.)  
(D x W x L)  
Part Number 99-9581